

# Safety Cabinets



For Ultra Clean Work Environment





## Safety Balance Work Station :

Safety Cabinets are designed for balance applications in the safety storage industry. Handling of harmful and toxic compounds are done safely thus reducing the risk of explosion. Safety Cabinets are high performance unit that are easy to use, clean and maintain. Particular attention is also paid to reliability and smooth running. Our commitment to customer needs has resulted in the most comprehensive line of safety storage products available.

### Features:

- ✓ Utility hole on the rear of hood.
- ✓ Acrylic enclosure.
- ✓ Back draft construction.
- ✓ Safe and comfortable work space.
- ✓ Easy to move with wheel caster.



### Technical Specifications & Ordering Information :

Model No.		LSBWS-01A	LSBWS-01B
Dimension (WxDxH) mm	Cabinet	1100 x 560 x 1345	900 x 560 x 1345
Body		Acrylic Sheet & SS 41 1.2t (epoxy powder coated)	
Table		SS 41 1.2t (epoxy powder coated) and aluminum bar	
Work Top		Marble stone 30t	
Catalog No.		25280101	25280102

### Accessories :

Model No.	TABLE		Filter Box FSBWS
	TSBWS (L)	TSBWS (S)	
Dimensions WxDxH (mm)	1100 x 500 x 720	900 x 560 x 720	500 x 500 x 850
Catalog No.	25280111	25080112	25280113



## Safety Trimming Work Station :

Safety Cabinets are specially designed for asepsis pertaining activities. Many critical applications in the medical, pharmaceutical and scientific fields demand an ultra clean work environment which is free from biological and particulate contamination. Cabinets are designed to provide a high degree of aseptic protection for products and apparatus in laboratory and production unit.

- Features:**
- ✓ Aseptic dispensing of media.
  - ✓ Cell culturing.
  - ✓ Preparation of medical devices and equipment.
  - ✓ Sterility testing.
  - ✓ Electronic valve for water and foot switch.
  - ✓ Easy to move with wheel caster.
  - ✓ Acrylic enclosure.
  - ✓ Back draft construction.



### Technical Specifications & Ordering Information :

Model No.		LSTWS-01A	LSTWS-01B
Dimension (WxDxH) mm	Cabinet	1100 x 550 x 1367	900 x 550 x 1367
Body		Acrylic Sheet & SUS 304 1.0t	
Table		SS 41 1.2t (epoxy powder coated)	
Work Top		Trespa 20t	
Cup Sink		SUS 304 1.0t	
Catalog No.		25280201	25280202

### Accessories :

Model No.	TABLE		Filter Box FSTWS
	TSTWS (L)	TSTWS (S)	
Dimensions WxDxH (mm)	1100 x 500 x 720	900 x 560 x 720	500 x 500 x 850
Catalog No.	25280221	25080222	25280223



4340 Christy St. Fremont, CA 94538 USA.  
Toll Free : (877) 620 9992  
Tel. : (925) 271 4322  
Fax : (925) 886 0400  
Email : [info@labnics.com](mailto:info@labnics.com)  
Website : [www.labnics.com](http://www.labnics.com)